

SISRI - Gastric Aspecific Reflex (SISRI. G. A. R.) mean intense digital pressure above the sternum midline or iliac crest or spleen projection of the spleen– (SISRI trigger points)

Latency time (Lt) in seconds	Latency time after preconditioning	MFR in seconds	fD & equilibria	SISRI (hyperfunction syndrome of lattice endothelial system)	Trigger points	Restano Manouvre (Apnea test for 5 seconds + Boxer test) Pre-clinical stress test	Diagnosis
Lt = 10	Lt = 20	3 < MFR < 4 normal MFR, associated activation, outcome +	fD ≥ 3 (ideal value fD=3.81) stange attractor	Spleen decongestion	Sternum Iliac crest Spleen	Lt = 9 seconds	Health
Lt < 10 [Lt = 3]	Lt < 20	MFR > 4 compromised MFR	1 < fD < 3 limit cycle to fix point	Complete SISRI	Sternum Iliac crest Spleen	Lt < 9 seconds	Infection by Positive GRAM (i.e. mumps) [Overt cancer]
Lt < 10	Lt < 20	MFR > 4 compromised MFR	1 < fD < 3 limit cycle to fix point	Intermediate SISRI	Sternum Iliac crest Spleen (less)	Lt < 9 seconds	Infection by Negative GRAM (Neisseria Meningitidis, Haemophilus Influenzae, Escherichia coli and Helicobacter pylori)
Lt < 10	Lt < 20	MFR > 4 compromised MFR	1 < fD < 3 limit cycle to fix point	Incomplete SISRI	Sternum Iliac crest	Lt < 9 seconds	Incomplete SISRI (i.e., virus infection)

Table 1. Legend: MFR (Microcirculatory Functional Reserve); EBD (Endoarteriolar Blocking Device); CAD (Coronary Artery Disease); fD (fractal Dimension); Lt (Latency time)